Limitations and restrictions:

The repetition of the indicated processes has a minimal effect on the galimplant® surgical drills. The useful life of these is normally determined by the wear and tear and damage caused due to repeated use.

The best way to assess the state of deterioration is by carrying out a thorough and functional check prior to use. Worn down cutting edges or signs of corrosion are indications that said drill is not suitable for use and it should be disposed of. The surgical drills should not be resharpened.

Preparation in the place of use

As soon as possible following their use, the galimplant® surgical drills should be transferred to the place where they are to be cleaned, in order avoid the drying-on of any blood residues, tissues, secretions or bone fragments. If it is likely that the cleaning process will be delayed, these should be covered with a moist cloth or stored in an airtight container filled with water in order to delay the drying-on of the organic material and to avoid the risk of contamination. To avoid contamination, the contaminated drills should not be put back into the surgical box. The used devices should be placed in a container that has been set up for the pre-cleaning which shall be carried out as per the following indications.

Preparation before cleaning

This stage must be carried out regardless as to whether the cleaning and disinfection processes are carried out manually or automatically.

- As soon as possible following their use, immerse the galimplant® surgical drills in a container that contains a cleansing solution with proteolytic enzymes, prepared in accordance with the manufacturer's specifications, respecting the temperature, concentration and immersion time. These enzyme solutions separate proteins and ensure that blood and other secretions do not dry on the devices.
- Make sure that all the drill's surfaces are covered by the solution and that the instruments do not touch each other.
- Manually brush with a nylon bristle brush (never metal, wool or steel brushes) until all of the visible dirt and residue has been eliminated, paying special attention to the zones which may not be reached by the brushing action, such as grooves or internal parts.
- Remove the drills from the container containing solution and rinse them thoroughly with plenty of water, preferably purified water, to ensure the reduction of organic material from mechanical drag.

Cleaning and disinfection

To clean and disinfect the drills either of the two following methods may be used:

Manual Method

- Place the galimplant® surgical drills in a perforated tray or basket.
- Prepare an ultrasonic bath and a detergent and disinfectant solution in accordance with the manufacturer's specifications of the product used. Nueva Galimplant has validated this process using the product neodisher® septo PreClean diluted at 1% to prepare the ultrasonic bath.
- Immerse the basket with the galimplant® surgical drills in the ultrasonic bath and keep them at a temperature of 20° C for 15 minutes.
- Remove the drills from the ultrasonic bath and rinse with plenty of purified water.
- If necessary, repeat the cleaning and disinfection step.
- Dry the devices using filtered compressed air and single use lint free cloth.

Automatic Method

Place the galimplant® surgical drills in the thermal washer-disinfector's tray or basket. Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.

- Run the cleaning, disinfecting and drying cycle.
- Nueva Galimplant has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:
 - Cold water soak for 4 minutes.
 - 55° C clean for 10 minutes with detergent.
 - Cold water neutralisation for 6 minutes.
 - Cold water mid-cycle rinse for 3 minutes.
 - 93° C thermal disinfection for 5 minutes.
 - 110° C dry for 25 minutes.

To guarantee optimum protection of the attachments we recommend that purified water is used throughout the whole process.

Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfector manufacturers.

Inspection, maintenance and testing

- Carefully inspect all of the galimplant® surgical drills in order to check if any damage has been caused to them (corrosion, loss of shape, breaking) or contamination. We recommend that the inspection is carried out in a well-lit environment or using a magnifying lens in order to improve visibility.
- If there are any traces of contamination on any of the devices the cleaning and disinfection process should be repeated.
- If there are any damage (deformation, failure of the drill, loss cutting capacity...), the drill should be discarded.

Packaging

Once they are clean, disinfected and completely dry, the drills must be packaged in a surgical box and then the box should be packaged in a sterilisation bag. The packaging of thermally sterilised medical devices must comply with the following requisites: - Be appropriate for steam sterilisation, that is, have a minimum thermic resistance of 137° C and be sufficiently permeable for steam.

- Protect the individual devices or surgical boxes sufficiently from mechanical damage.

- Comply with the EN ISO 11607 normative.

Once the devices have been packaged a label should be added to each sterilisation bag which contains data such as the sterilisation date and it should include the reference, batch and the number of devices if the bag does not have a transparent window which allows the user to see what is inside.

Sterilisation

- Load the autoclave with the conditioned galimplant® surgical drills in the sterilisation bags. Do not overload the autoclave and always follow the equipment manufacturer's instructions.
- Just like in the other steps of these instructions, we recommend the use of purified water to guarantee optimum quality.
- Nueva Galimplant recommends using an autoclave which complies with the requisites of the EN ISO 17665 normative, in addition to being validated and having adequate maintenance. The validation of the sterilisation was carried out in a steam autoclave using a 134° C pre-vacuum cycle for 6 minutes.
- To control the sterilisation process, the use of chemical indicators is recommended, such as, for example, the verification strips which in the sterilisation bags come equipped with and which change colour when the sterilisation process has finished.

Note: it is important to always follow the autoclave manufacturer's instructions and keep maintenance up to date.

WARNINGS

Use the surgical drills only for the insertion of original galimplant® dental implants.

The practitioner must provided the patient with the information about side effects, infections and complications in relation to galimplant® surgical drills.

The clinical practitioner must verify that the materials and techniques to be used during procedure are appropriate and compatible with the galimplant® surgical drills and must follow manufacturer's instructions for use of these materials.

Failure to observe the procedures set out in these instructions could harm the patient and give rise to one or all of the following complications:

- Aspiration or swallowing of the components.
- Infection.
- Damage to the implant or instruments.
- Implant failure.

In the event of a serious incident related to the use or result of use of this product, inform Nueva Galimplant (www.galimplant.com) and your National Competent Authority.

Procedures for the disposal of the product

If used, dispose the surgical drills as biological waste, according to the local regulations.

Validity

These instructions for use will be valid as of May 2022.

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