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IMPORTANT! Read the information contained in these operating instructions carefully before using the product.

NUEVA GALIMPLANT complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical products. The galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the galimplant® system is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U.

The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations.

Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their obligation to personally make sure that the product is appropriate for the purposes and procedures planned.

Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about its products through the platform: galimplant.training

Precautions and warnings

In general terms, our products must be held to avoid their aspiration during intraoral use or insertion in the implant.

Galimplant® multi-position angled abutments have not been sterilised before delivery, so they must be cleaned, disinfected and sterilised before use and it is essential that they are handled in an aseptic manner.

Galimplant® multi-position angled abutments are of single use. The reusing of these types of products could lead to a deterioration of their characteristics which involves risk of tissue infection, surgical or prosthodontic failure and/or damage to the patient's health.

The following descriptions are insufficient for the immediate application of any abutment from the galimplant® system, this system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of galimplant® products, in order to use them safely and correctly, following the indications of these instructions for use.

Description

Galimplant® multi-position angled abutments are made from grade 5 titanium and can be used in galimplant® implants with different connection types (external and internal), platforms, diameters and lengths.

See in the PRODUCTS section the implant models compatible with each of the configurations presented below.

Our multi-position angled abutments are available in the following configurations:

- Rotational Internal Connection:
 - Multi-position abutment angled at 15° with shoulder heights of 2 and 3 mm and pink colour anodized.
 - Multi-position abutment angled at 30° with shoulder height of 3 mm and pink colour anodized.
 - Multi-position abutment angled at 45° with shoulder height of 5 mm and pink colour anodized.
 - Multi-position abutment angled at 50° with shoulder height of 5 mm and pink colour anodized.
 - Multi-position abutment angled at 60° with shoulder height of 5 mm and pink colour anodized.
- Rotational External Connection:
 - Multi-position abutment angled at 15° with shoulder heights of 2 and 3 mm and pink colour anodized.
 - Multi-position abutment angled at 30° with shoulder height 3 and 4 mm and pink colour anodized.

The shoulder height will depend on the soft tissue height and the angulation will depend on the degree of the implant's inclination.

They come with:

- 1 transporter
 - 2 screws blue colour anodized (one clinical, packed inside a capsule and another prosthetic)
 - 1 transeptal transfer attachment
 - 1 screw for the transeptal transfer attachment
 - 1 transeptal closing abutment
 - Recommended torque:
 - 30 Ncm screw to fix the abutment to the implant
 - 10 Ncm prosthetic screw to fix the abutment to the prosthesis
- Do not reuse. Single use.

Intended use

Our multi-position angled abutments with direct implant connection are designed for multi-unit screw-retained restorations, with the aim of restoring the patient's mastication function and phonetics.

Indications

Our multi-position angled abutments are recommended in the following cases:

- Fixed multi-unit screw-retained restorations.
- Restorations of the posterior sector.
- Correction of implants with large non-parallelisms.

Patient Groups and Clinical Benefits

The galimplant® multi-position angled abutments with direct implant connection are intended for use in totally or partially edentulous people or with poor prognosis or non-restorable teeth that request treatment and that meet the criteria for psychological, medical, anatomical, surgical, prosthodontic and maintenance selection specified in these instructions for use.

The clinical benefit of treatment is the replacement of the patient's missing teeth.

Contraindications

The use of multi-position angled abutments is not recommended for patients:

- who do not meet the necessary medical requirements for dental implant placement surgery.
- who are allergic or hypersensitive to commercially pure titanium or Ti-6Al-4V titanium alloy (titanium, aluminium and vanadium).
- Restoration of the anterior sector.
- When the definitive orifice of the screw entry is in an area of elevated aesthetic compromise.

Side effects, interactions and complications

The side effects, interactions and complications are those derived from implant placement procedures.

Possible complications following prosthetic restoration include:

- Complications: Loosening/fracture of the abutment and/or screw, mechanical problems associated with the prosthesis, loss of prosthesis retention, loss of the hole sealing of the screw-retained prosthesis and microgap in the implant-abutment connection.

If one of these disorders/complications appears, the patient should consult their specialist doctor who will give them the corresponding instructions.

Storage and transport

The multi-position angled abutments are delivered packed in peel-pack bags inside a closed cardboard blister on which the label is placed and includes all product information, including the date of manufacture, which can be deduced from the batch number, the first four digits corresponding to the year and the next two to the month of manufacture (YYYYMMXXX). This packaging contains and protects the product during storage and transport and does not have a specified shelf-life.

The multi-position angled abutments must be stored in a dry place in the original packaging at room temperature and not exposed to direct sunlight.

When removing the product from its packaging one must verify that it is not damaged and that the product fits the description on the label.

Cleaning/Sterilisation/Conservation

Caution! All our multi-position angled abutments are sold in a NON STERILE state.

Before use, they must be cleaned, disinfected and sterilised, observing the following procedures validated by Nueva Galimplant SLU. Said procedures should be carried out before intraoral use. The repetition of the procedures described in this section do not alter the characteristics of these devices. The non-observance of these instructions may cause cross infections.

The use of neutral or slightly alkaline cleaning agents is recommended. It is essential to verify that the cleaning product used is compatible with the medical devices are going to

processed and respect the specifications of use (temperature, concentration and time) indicated by the manufacturer.

Preparation before cleaning

- Immerse all the attachments in a container that contains a cleansing solution with proteolytic enzymes, prepared in accordance with the manufacturer's specifications, respecting the temperature, concentration and immersion time.
- Make sure that all the attachments' surfaces are covered by the solution and that the attachments do not touch each other.
- Manually brush with a nylon bristle brush (never metal, wool or steel brushes) until eliminating all the visible dirt and residue, paying special attention to the zones which may not be reached by the brushing action, such as grooves or internal parts.
- Remove the attachments from the container containing solution and rinse them thoroughly with plenty of water, preferably purified water, to ensure the reduction of organic material from mechanical drag.

Cleaning and disinfection

To clean and disinfect the attachments either of the two following methods may be used:

Manual Method

- Place the attachments in a perforated tray or basket.
- Prepare an ultrasonic bath and a detergent and disinfectant solution in accordance with the manufacturer's specifications of the product used. *Nueva Galimplant S.L.U. has validated this process using the product neodisher® septo PreClean diluted at 1% to prepare the ultrasonic bath.*
- Immerse the basket with the attachments in the ultrasonic bath and keep them at a temperature of 20° C.
- Remove the attachments from the ultrasonic bath and rinse with plenty of purified water.
- If necessary, repeat the cleaning and disinfection step.

Automatic Method

- Place the attachments in the thermal washer-disinfector's tray or basket. Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.
- Run the cleaning, disinfecting and drying cycle. *Nueva Galimplant S.L.U. has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:*
 - Cold water soak for 4 minutes.
 - 55° C clean for 10 minutes with detergent.
 - Cold water neutralisation for 6 minutes.
 - Cold water mid-cycle rinse for 3 minutes.
 - 93° C thermal disinfection for 5 minutes.
 - 110° C dry for 25 minutes.

To guarantee optimum protection of the attachments we recommend that purified water is used throughout the whole process.

Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfector manufacturers.

Packaging

Once they are clean, disinfected and completely dry, the attachments must be packaged in a sterilisation bag.

The packaging of thermally sterilised medical devices must comply with the following requisites:

- Be appropriate for steam sterilisation, that is, have a minimum thermic resistance of 137° C and be sufficiently permeable for steam.
- Protect the devices sufficiently from mechanical damage.
- Comply with the EN ISO 11607 normative.

Sterilisation

- Load the autoclave with the conditioned attachments in the sterilisation bags. Do not overload the autoclave and always follow the equipment manufacturer's instructions.
- Just like in the other steps of these instructions, we recommend the use of purified water to guarantee optimum quality.
- Nueva Galimplant S.L.U. recommends using an autoclave which complies with the requisites of the EN ISO 17665 normative, in addition to being validated and having adequate maintenance. The validation of the sterilisation was carried out in a steam autoclave using a 134° C pre-vacuum cycle for 6 minutes.
- To control the sterilisation process, the use of chemical indicators is recommended, such as, for example, the verification strips which the sterilisation bags come equipped with and which change colour when the sterilisation process has finished.

Note: it is important to always follow the autoclave manufacturer's instructions and keep maintenance up to date.

Caution! The re-use and/or re-sterilisation of prosthetic abutments can cause risk of infection and/or cross contamination.

Documentation, traceability and implant card.

Nueva Galimplant S.L.U. recommends keeping an exhaustive clinical, radiological and photographic documentation. All products can be identified through their reference and batch code and / or UDI code. To maintain the mandatory traceability of these products, two adhesive labels are included, one for the medical record and the other for the patient, the latter must be attached to the implant card supplied with the products and which must be delivered to the patient fully completed once the surgery is over.

WARNINGS

Only original prosthetic abutments should be placed in the galimplant® implants. The prosthetic and clinical practitioner must verify that the materials and techniques to be used during prosthetic and clinical procedures are appropriate and compatible with the galimplant® multi-position angled abutments and must follow the manufacturer's instructions for use of these materials.

Before the prosthetic part, the physician should carry out an exhaustive patient study, in order to determine a correct planning of the case.

For the placement of the multi-position angled abutments, use only original wrenches of the brand galimplant®. Consult the compatible product models in the catalogue available at www.galimplant.com.

There is a risk of aspiration or ingestion of the products when they are used, for this reason, they should only be used by professionals.

The prosthetic abutments must solely and exclusively be placed during the implants osteointegration or if the implants have a great primary stability.

Failure to observe the procedures set out in these instructions could harm the patient and give rise to one or all of the following complications:

- Aspiration or swallowing of the components.
- Infection.
- Damage to the implant, abutment, components or instruments.
- Implant failure.

To guarantee a good long term treatment result, the clinician must inform the patient about appropriate home oral hygiene and the importance of attending periodic clinical and radiographic consultations to perform a control and preventive diagnosis of any side effect and/or complication.

In the event of a serious incident related to the use or result of use of the product, inform Nueva Galimplant (www.galimplant.com) and your National Competent Authority.

The practitioner must provided the patient with the information about side effects, interactions and complications in relation to galimplant® multi-position angled abutments.

Procedures for the disposal of the product

Once they have been extracted from the oral cavity, the prosthetic components must be disposed of as biological waste, in accordance with current local laws.

Summary of Safety and Clinical Performance

The Summary of Safety and Clinical Performance (SSCP) for the product is available in the European database EUDAMED. The document can be obtained from the website: <https://ec.europa.eu/tools/eudamed>.

Validity

These instructions for use will be valid as of May 2022.

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