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English

IMPORTANT! Read the information included in these instructions for use before using the product

NUEVA GALIMPLANT complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical products. The Galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the Galimplant® System is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions.

This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U.

The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations.

Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from his obligation to personally make sure that the product is appropriate for the purposes and procedures planned.

Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about its products through the platform: galimplant.training.

Precautions and warnings

The galimplant® dental implants are delivered sterile, they are for single use and must not be used after their expiration date.

Warning! Re-sterilisation and/or reuse of implants may result in the risk of crossed contamination.

The following descriptions are insufficient for the immediate application of any Galimplant® system implants. This system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of the product in order to use them safely and correctly following the indications of these instructions for use.

Description

The internal connection dental implants with implant transfer for closed tray technique are supplied fitted in a plastic implant transfer attached to the bone, which can be used to impress taking with the closed tray technique. They can be supplied with and without closing abutment. Refer to the instructions for use without closing abutment are those ending in S (see table attached at the end of this instruction).

The endosseous implant is manufactured using commercially pure grade IV titanium, in accordance with active ISO 5832-2 and presents a Nanoblastium® surface, the result of a technique which combines a jet of sand and hydrofluoric acid etching. The closing abutment and implant transfer are manufactured using grade V titanium (Ti6Al4V) in accordance with ISO 5832-3.

The internal connection dental implants with implant transfer for closed tray technique are available with internal hexagonal connection, manufactured in turn, in different models, diameters and lengths.

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The selection of the model, diameter and length will depend on the height and width of the bone, the bone density and the surgical technique which is to be used. Check out the implant models available in the PRODUCTS section.

Compatibility

The ICIC and IPXC implant models are compatible with galimplant® abutments and screws with internal hexagonal connections.

Galimplant® surgical drill system must be used to prepare the implant bed, following the indicated drill system and only the original galimplant® wrenches should be used for the handling and placement of implants, as stated in the instructions for use.

Intended use

The galimplant® implant system is designed for use in oral implantology in order to re-establish lost functions.

Indications

The dental implants are indicated for use in partial or complete restoration for patients who are partially or completely edentulous and they may be used in all of the indications for oral endosseous maxillary and mandibular implants.

Once the dental implants have been osseointegrated serve as a support for different fixed or removable prosthetic solutions that allow for the restoration of the patient's mastication and phonetics function.

Patient Groups and Clinical Benefits

The galimplant® implant system is intended for use in totally or partially edentulous people or with poor prognosis or non-restorable teeth that request treatment and that meet the criteria for psychological, medical, anatomical, surgical, prosthodontic and maintenance selection specified in these instructions for use. The clinical benefit of treatment is the replacement of the patient's missing teeth.

Contraindications

The applicable contraindications are those related to surgical procedures for the fitting of implants.

-Absolute contraindications: Severe internal medical problems, metabolic bone disorders, uncontrolled coagulopathy, un-cooperative or un-motivated patients, alcoholism or drug addiction, psychosis, history of treatment resistant disorder, xerostomia, weakened immune system, illnesses treated with regular usage of steroids, allergy or hypersensitivity to titanium, uncontrollable endocrine illnesses.

- Relative contraindications, in these specific circumstances in which the risk may not be considered as high as the benefit obtained after implant placement, the clinician must assess these conditions: Previously irradiated bone, diabetes, anticoagulation therapy/hemorrhagic diathesis, bruxism, parafunctional habits, unfavourable bone anatomy, smoking, uncontrolled periodontitis, temporomandibular joint disorder, maxillary disease and altered oral mucous membrane susceptible to treatment, controllable endocrine illnesses, pregnancy, lactation, children with unfinished bone growth, insufficient oral hygiene, insufficient bone available or inadequate bone quality, localized radicular remains.

Side effects, interactions and complications

The side effects, interactions and complications are those derived from implant placement procedures. Possible complications following dental implant placement procedures include:

- Temporary disorders: Pain, inflammation, difficulty speaking, gingivitis, edema, swelling, difficulty opening the mouth, bad breath, labial laceration, haemorrhage, contusion, hematoma, temporal paraesthesia.
- Longer-lasting problems: Chronic pain associated with the dental implant, permanent paraesthesia, dyesthesia, bone loss from the maxillary crest, localised or systemic infection, oral-antral or oral-nasal fistulas, periimplantitis, implant failure, irreversible damage to neighbouring teeth, implant, maxillary or prosthetic fractures, aesthetic problems, nerve damage, exfoliation, hyperplasia, periapical infection, mucositis, marginal bone loss, sinusitis, perforation of the anatomical areas.

In case in which any of these disorders appear, the patient should report them to the specialist who will give the corresponding indications.

Packaging and sterility

The dental implants have been manufactured and packaged in a controlled environment prior to their sterilisation by gamma rays.

The metallic implant transfer that holds the endosseous implant and the closing abutment are presented assembled on a plastic support. The aforementioned plastic support is presented inside a container that guarantees sterility by means of an opérculo.

Both the external and the internal packaging bears a label which indicates the following concepts: implant model, size, batch and expiration date, in addition the external label seals the packaging guaranteeing that the content cannot be manipulated.

Nueva Galimplant S.L.U. decline all responsibility in the case of re-sterilizing implants, irrespective of the person performing the sterilization or the method used. A previously used or non-sterile implant should not be implanted under any circumstance. Nueva Galimplant S.L.U. cannot accept the return of an implant if the packaging is damaged.

Storage and transport

The implants should be maintained in their original packaging, in a dry environment, at room temperature and away from direct sunlight. These same conditions are also required during transportation and in addition the product should not be knocked or crushed.

The intact sterile packaging protects the implant from any external influences during storage and transportation, and if it is stored in adequate conditions it guarantees the sterility until the expiration date. The implant should not be used following the expiration date (printed on the label).

When taking the product out, ensure that there is no damage and that the product is as described on the label.

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