IMPLANTE CIGOMÁTICO GALIMPLANT GALIMPLANT ZYGOMATIC IMPLANT IMPLANTE ZIGOMÁTICO GALIMPLANT

Español

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English

IMPORTANT! Read the information included in these instructions for use before using the product

NUEVA GALIMPLANT complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical products. The Galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the Galimplant® System is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-

recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U. The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations.

Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their obligation to personally make sure that the product is appropriate for the purposes and procedures planned.

Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about it products through the platform: galimplant.training.

Precautions and warnings

The galimplant® dental implants are delivered sterile, they are for single use and must not be used after their expiration date.

Warning! Re-sterilisation and/or reuse of implants may result in the risk of crossed contamination.

The following descriptions are insufficient for the immediate application of any Galimplant® system implants. This system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of the product in order to use them safely and correctly following the indications of these instructions for use.

Description

The zygomatic implants are supplied fitted in a titanium implant transfer connected to the implant connection. This implant transfer can be used as an impression coping for the closed tray technique. These can be supplied with and without closing abutment. References without a closing abutment are those ending in S (see table attached at the end of this instruction).

The endosseous implant is manufactured using commercially pure grade IV titanium, in accordance with directive ISO 5832-2 and it has a Nanoblast^{plus} surface, the result of a technique which combines a jet of sand and a triple acid etching. The closing abutment, implant transfer and implant transfer screw are manufactured using grade V titanium (Ti6Al4V) in accordance with ISO 5832-3.

The zygomatic implants are available with two types of design, manufactured in turn, in different diameters and lengths.

- **Smooth Zygomatic Implant (ICMS)**: designed with hexagonal internal connection, it has a threaded apex with Nanoblast^{plus}® treatment and machined smooth base with a 4.2mm platform and is available in lengths of 30 to 60mm at 2.5mm intervals.
- **Treated Zygomatic Implant (ICMT)**: designed with hexagonal conical internal connection, it has an apex thread, a smooth threaded neck and a serrated surface with Nanoblast^{plus}® treatment, it has a 4.3mm platform and is available in lengths of 20 to 60mm at 2.5mm intervals.

The selection of design and length will depend on the anatomy, the quality and quantity of bone and the space available in each clinical case. Check the available implant models in the PRODUCTS section.

Compatibility

The zygomatic implants are compatible with galimplant® fitted multi-position abutments with internal hexagonal connection, available with angles 15°, 30°, 45°, 50° and 60°. For the preparation of the impression of the galimplant® surgical drill system must be used, respecting the indicated drilling sequence and the original galimplant® techniques indicated in the procedure for use that should be used for the handling and placement of implants.

The zygomatic implants are endosseous implants to integrate into the zygomatic bone and are designed for surgical placement in the bone of the upper jaw to support prosthetic devices, such as artificial teeth, in order to restore the mastication function of the phonetics of the patient.

Indications

The zygomatic implants are indicated in cases with severe maxillary atrophy for the restoration of totally or partially edentulous patients, and can be used in combination with other implants of the Galimplant® Sterile Dental Implant System family.

Patient Groups and Clinical Benefits

The galimplant® zygomatic system is intended for use in totally or partially edentulous people or with poor prognosis or non-restorable teeth that request treatment and that meet the criteria for psychological, medical, anatomical, surgical, prosthodontic and maintenance selection specified in these instructions for use.

The clinical benefit of treatment is the replacement of the patient's missing teeth.

Contraindications

The applicable contraindications are those related to surgical procedures for the fitting of implants. -Absolute contraindications: Severe internal medical problems, metabolic bone disorders, uncontrolled osteoporosis, un-cooperative or un-motivated patients, alcoholism or drug addiction, psychosis, history of treatment resistant disorder, dry mouth, weakened immune system, illnesses treated with regular use of steroids, allergy or hypersensitivity to titanium, uncontrollable endocrine illnesses.

-Relative contraindications, in these specific circumstances in which the risk may not be considered as high as the benefit obtained after implant placement, the clinician must assess these conditions: previously irradiated bone, diabetes, anticoagulation therapy/hemorrhagic diathesis, bruxism, parafunctional habits, unfavourable bone anatomy, smoking, uncontrolled periodontitis, temporomandibular joint disorder, maxillary disease and altered oral mucous membrane susceptible to treatment, controllable endocrine illnesses, pregnancy, lactation, children with unfinished bone growth, insufficient oral hygiene, insufficient bone available or inadequate bone quality, localized radicular remains.

In addition to these contraindications related to the surgical procedures for implant placement described above, the zygomatic implants are contraindicated in single rehabilitations and in patients with any pathology of the zygomatic bone and sinus.

Side effects, interactions and complications

The side effects, interactions and complications are those derived from implant placement procedures. Possible complications following dental implant placement procedures include:

- Temporary disorders: Pain, inflammation, difficulty speaking, gingivitis, edema, swelling, difficulty opening the mouth, bad breath, labial laceration, haemorrhage, contusion, hematoma, temporal paraesthesia.
- Longer-lasting problems: Chronic pain associated with the dental implant, permanent paraesthesia, dyesthesia, bone loss from the maxillary crest, localised or systemic infection, oral-antral or oral-nasal fistulas, periimplantitis, implant failure, irreversible damage to neighbouring teeth, implant, maxillary or prosthetic fractures, aesthetic problems, nerve damage, exfoliation, hyperplasia, periapical infection, mucositis, marginal bone loss, sinusitis, ocular injury, paraesthesia of the infraorbital and facial zygomatic nerves, perforation of the anatomical areas.

In case in which any of these disorders appear, the patient should report them to the specialist who will give the corresponding indications.

Packaging and sterility

The zygomatic implants have been manufactured and packaged in a controlled environment prior to their sterilisation by gamma rays.

The metallic implant transfer that holds the endosseous implant and the closing abutment are presented assembled on a plastic support. The aforementioned plastic support is presented inside a container that guarantees sterility by means of an opérculo.

Both the external and the internal packaging bears a label which indicates the following concepts: implant model, size, batch and expiration date, in addition the external label seals the packaging guaranteeing that the content cannot be manipulated.

Nueva Galimplant S.L.U. decline all responsibility in the case of re-sterilizing implants, irrespective of the person performing the sterilization or the method used. A previously used or non-sterile implant should not be implanted under any circumstance. Nueva Galimplant S.L.U. cannot accept the return of an implant if the packaging is damaged.

Storage and transport

The implants should be maintained in their original packaging, in a dry environment, at room temperature and away from direct sunlight. These same conditions are also required during transportation and in addition the product should not be knocked or crushed.

The intact sterile packaging protects the implant from any external influences during storage and transportation, and if it is stored in adequate conditions it guarantees the sterility until the expiration date. The implant should not be used following the expiration date (printed on the label).

When taking the product out, ensure that there is no damage and that the product is as described on the label.

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Planing

The candidate patient to be a carrier of zygomatic implants must not present any medical / surgical contraindication that discourages their use.

Before surgery, an exhaustive patient evaluation and study should be carried out (assessment, study models, panoramic and periapical x-rays and C.A.T. scan) in order to determine the anatomical characteristics, occlusion conditions, periodontal condition and to determine the bone suitability in order to draw up the treatment plan. The physician is exclusively authorised to plan this treatment. Faults in the treatment plan may be the cause of complications related to zygomatic implants: either prosthetic, mechanical, biological and may even lead to the loss of the implant(s).

Surgical Procedure

Surgery for the placement of the galimplant® zygomatic implants can be performed under the following anesthetic techniques: local anesthesia, intravenous sedation or general anesthesia. The combination of different anesthetic procedures can be used.

To begin exposure of the lateral maxillary wall, a full thickness mucoperiosteal flap is reflected following a crestal incision, which displaces at least 3 mm of keratinized gingiva towards vestibular. The crestal incision will be accompanied by a mesial and / or distal vertical discharge, as appropriate, to allow exposure of a tension-free operative field and with adequate exposure of the area to be treated.

- It starts in the alveolar portion, either with round or extra-long lanceolate drill at speeds between 600-1200 rpm, performing the drilling with abundant and constant irrigation of sterile saline.
- Once the proposed route has been established, the implant bed will be shaped according to a drilling sequence with gradual increments of diameters (drills of diameters 1, 2 and 3). For bone densities in the body of the malar bone type D4, the use of the first drill - number 1 - will suffice; in cases of densities type D2-D3, it will be necessary to pass drills number 1 and 2; and for the bone grades Type D1, the drills 1 and 2 will be used up to the final distance of working, and the drill number 3 at 5 mm less than the final distance (easily recognizable by the marks found on them). All the drills and components are marked to prepare the implant bed at the correct depth and achieve a safe and predictable position.
- The shape for housing the implant body can be made with the side cutting drill. Depending on the anatomy of the patient, the implant body can be placed inside or outside the maxillary sinus. Ensure correct angulation and prevent the drills from wobbling, as this may unintentionally widen the preparation bed. During the entire drilling phase, profuse irrigation of both the alveolar portion and the body of the zygomatic bone and the lateral wall of the sinus is recommended.

• Implant length measurement

Use the depth gauges to determine the length of the zygomatic implant to be placed. These gauges reflect, at 5 mm intervals, the different measurements available.

Note: It must be taken into consideration that implant measurements do not consider the prosthetic abutment. For 45°, 50° and 60° angulations, the length from the head of the implant to the center of the abutment corresponds to 5 mm.

• Placement of zygomatic implants

All zygomatic implants are supplied with an implant transfer retained by friction. The placement of the implant, after the preparation of the implant bed; it can be done manually, using the appropriate wrenches (long or short driver) or using the manual screwdriver of the surgical box. To take the implant out of the container open the external packaging, unscrew the plastic top and remove the operculum.

The implant may be taken out of the package in the following ways:

- Connect the short or long driver to the metal implant transfer.
- Placing the manual screwdriver tip directly on the implant transfer.

Once connected remove the implant from the plastic holder by applying light pressure in an upwards direction at the aperture where it is placed and hold it with any of the assembly techniques indicated above to bring it to the mouth and insert it into the implant bed. Once the implant has been inserted, position it according to the desired depth, not exceeding an insertion torque greater than 50 Ncm.

Once the implant is placed in position verify that the stability is adequate (but does not exceed 50 Ncm), that no more than 2 mm emerge from the side wall of the zygomatic bone, that it does not invade any noble structure and that the placement of the abutment prosthetic does not interfere with any other implant / abutment or dental unit.

If the case requires the placement of conventional implants, they are placed according to the conventional protocol.

• Selection of intermediate abutment

Nueva Galimplant recommends the selection of the transepithelial abutment at the time of surgery, since it will only be at this time when can be corrected both the depth of insertion and the prosthetic emergence perpendicular to the occlusal plane and parallel to the rest of the implant/s if any. To facilitate the selection of the same, there is a set of guides that determine the required angulation. The abutments (multi-position) available, present the following angulations: 0°, 15°, 30°, 45°, 50° and 60°. All the abutments should be torqued to 50 Ncm.

For immediate loading, all implants must have an insertion torque greater than 35 Ncm., must have a placement that allows a proper distribution of forces, the abutments must be properly inserted and tightened, the provisional prosthesis should be sufficiently rigid, the occlusion will be conveniently established, as well as a strict revision plan should be established.

If the immediate loading can not be carried out, it is recommended to install the close abutment (5 Ncm), which is supplied with the implant anchored in the lower part of the plastic support. For the two-phase protocol, relieve the denture on the implants.

• Flap closure

Before proceeding to the adequate closure of the treated area, a profuse washing of the operative field is recommended, it will be conveniently checked that there are no bleeding areas, remains of surgical and / or prosthetic material, as well as material used during the surgery.

After implant placement, the surgeon's evaluation of the amount of bone and primary stability will determine if the implant loading can be performed.

Documentation, traceability and implant card.

Nueva Galimplant S.L.U. recommends keeping an exhaustive clinical, radiological and photographic documentation. All products can be identified through their reference and batch code and / or UDI code. To maintain the mandatory traceability of these products, two adhesive labels are included, one for the medical record and the other for the patient, the latter must be attached to the implant card supplied with the products and which must be delivered to the patient fully completed once the surgery is over.

WARNINGS

This product requires specific knowledge of the zygomatic implant placement technique. The physician should provide the patient with information regarding the contraindications, side effects, interactions and complications associated with the galimplant® implants.

The patient should avoid activities which require a high level of physical effort immediately after the insertion of the dental implant.

When using our products, the necessary measures should be taken to avoid inhalation.

Before placing the implant ensure that there is no damage has been caused to both containers. A damaged container may mean that the implant is now not sterile. When removing the implant from the sterile container follow the corresponding aseptic principles. The container should only be opened when the surgical intervention is to be carried out.

The health professional must ensure that the materials and techniques that are to be used during the clinical and prosthetic procedures are appropriate and compatible with the zygomatic implants and the instructions for use provided by the manufacturer of these materials should be followed.

For the preparation of the implant bed, and the placement and restoration of the zygomatic implants only original galimplant® products must be used.

Failure to follow the procedures summarised in these instructions may cause damage to the patient, may to one or more of the following complications.

- Inhalation or swallowing of the components.
- Infection
- Damage to the implant, abutment, components or instruments.

In the event of a serious incident related to the use or result of use of the product, inform Nueva Galimplant (www.galimplant.com) and your National Competent Authority.

Procedures for the disposal of the product

Once extracted from the oral cavity the prosthetic components must be disposed of as biological waste in accordance with the local regulations currently in force.

Summary of Safety and Clinical Performance

The Summary of Safety and Clinical Performance (SSCP) for the product is available on the EUDAMED database EUDAMED. The document can be obtained from the web page <https://ec.europa.eu/tools/eudamed>.

Validity

These instructions for use will be valid as of September 2023.

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