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English

IMPORTANT! Read the information contained in these operating instructions carefully before using the product.

NUEVA GALIMPLANT complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical products. The galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the galimplant® system is not under the manufacturer’s control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U.

The user has sole responsibility for the procedure and use of the product. It is the practitioner’s responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient’s individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations.

Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their obligation to personally make sure that the product is appropriate for the purposes and procedures planned. Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about it products through the platform: galimplant.training.

Precautions and warnings

In general terms, our products must be held to avoid their aspiration during intraoral use or insertion in the implant.

The galimplant® impressions coping have not been sterilised before delivery, so they must be cleaned, disinfected and sterilised before use and it is essential that they are handled in an aseptic manner.

The galimplant® impressions coping are of single use. The reusing of these types of products could lead to a deterioration of their characteristics which involves risk of tissue infection, surgical or prosthodontic failure and/or damage to the patient’s health.

The following descriptions are insufficient for the immediate application of any attachment from the galimplant® system, this system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of galimplant® products, in order to use them safely and correctly, following the indications of these instructions for use.

Description

The galimplant® impression copings are made from grade 5 titanium and can be used in galimplant® implants with different connection types (internal and external), platforms, diameters and lengths.

They have two marks with a distance of 1mm between them to help determine the available gingival height.

See in the PRODUCTS section the implant models compatible with each of the configurations presented below.

The impression coping screws are available in two lengths: 13.5 mm and 20 mm. Manufactured in grade V titanium.

The impression coping is available in the following configurations:

- Hexagonal IC impression coping supplied with short screw for impression coping (metrics 1,6mm).
- Hexagonal IC impression coping supplied with long screw for impression coping (metrics 1,6mm)
- Universal EC impression coping supplied with short screw for impression coping (metrics 2mm)
- Universal EC impression coping supplied with long screw for impression coping (metrics 2mm)
- Impression coping for rotational multi-position abutment supplied with short screw for impressions coping for multi-position (metrics 1,4mm)
- Non-rotational impression coping for multi-position abutments supplied with short screw for impressions coping for multi-position (metrics 1,4mm)
- Non-rotational impression coping for multi-position straight aesthetic supplied with short screw for impressions coping (metrics 2mm)
- Rotational impression coping for multi-position straight aesthetic supplied with short screw for impressions coping (metrics 2mm)
- Non-rotational impression coping for multi-position straight aesthetic supplied with long screw for impressions coping (metrics 2mm)
- Rotational impression coping for multi-position straight aesthetic supplied with long screw for impressions coping (metrics 2mm)
- Hexagonal IC impression coping for closed tray (implant transfer)

Do not re-use. Single use.

Intended use

The galimplant® impressions coping are provided solely and exclusively to transfer the position of the implant in the mouth for the study model on which the rehabilitation will be made.

Patient Groups and Clinical Benefits

The galimplant® impression copings are intended for use in patients subject to dental implant treatment.

The clinical benefit of treatment is the replacement of the patient’s missing teeth.

Contraindications

The use of galimplant® impressions copings are not recommended for patients:

- who do not meet the necessary medical requirements for dental implant surgery.
- who are allergic or hypersensitive to commercially pure titanium or Ti-6Al-4V titanium alloy (titanium, aluminium and vanadium).

Side effects, interactions and complications

The side effects, interactions and complications are those derived from implant placement procedures.

Storage and transport

The galimplant® impressions copings are delivered packed in peel-pack bags inside a closed cardboard blister on which the label is placed and includes all product information.

This packaging contains and protects the product during storage and transport and does not have a specified shelf-life.

The galimplant® impression copings are made of with stable materials under normal environmental conditions, but it is recommended to store them in the original packaging and not expose them to extreme conditions (temperatures, humidity, shock, etc.) that may affect the packaging and labeling.

When removing the product from its packaging one must verify that it is not damaged and that the product fits the description on the label.

Cleaning/Sterilisation/Conservation

The galimplant® impression copings are sold in a non-sterile state. In accordance with the EN ISO 17664 Normative, Nueva Galimplant will provide all of the information necessary to process their health products labelled as “Non sterile” and which need to be cleaned, disinfected and sterilized prior to use, both for the first use and for all subsequent usages. All models of the galimplant® impression coping are included in said health products.

Warnings:

- The new “non-sterile” products should be processed in the same way as the used products.
- We recommend carrying out this process as soon as possible following use.
- For their own safety, staff should use adequate protective clothing and equipment at all times: gloves, mask, protective glasses, etc;
- Nueva Galimplant recommends the use of aldehyde free cleaner and disinfectants with proven effectiveness, in order to avoid the fixation of blood proteins (included on the VAH/DGHH list or which have the CE Marking).
- Nueva Galimplant cannot control the clinical manipulation procedures, the cleaning methods or any other user level conditions, therefore they will not be held liable for the sterilization of the products even when the recommended guidelines contained in this instruction manual have been followed.

Limitations and restrictions

The repetition of the indicated processes has a minimal effect on the impression coping. The useful life of these is normally determined by the wear and tear and damage caused due to repeated use.

The best way to assess the state of deterioration is by carrying out a thorough and functional check prior to use.

Preparation in the place of use

As soon as possible following their use, the impression coping should be transferred to the place where they are to be cleaned, in order avoid the drying-on of any blood residues, tissues, secretions or bone fragments. If it is likely that the cleaning process will be delayed, these should be covered with a moist cloth or stored in an airtight container filled with water in order to delay the drying-on of the organic material and to avoid the risk of contamination.

Preparation before cleaning

- Immerse all the attachments in a container that contains a cleansing solution with proteolytic enzymes, prepared in accordance with the manufacturer’s specifications, respecting the temperature, concentration and immersion time. These enzyme solutions separate proteins and ensure that blood and other secretions do not dry on the devices.
- Make sure that all the attachment’s surface are covered by the solution and that the attachments do not touch each other.
- Manually brush with a nylon bristle brush (never metal, wool or steel brushes) until eliminating all the visible dirt and residue, paying special attention to the zones which may not be reached by the brushing action, such as grooves or internal parts.
- Remove the attachments from the container containing solution and rinse them thoroughly with plenty of water, preferably purified water, to ensure the reduction of organic material from mechanical drag.

Cleaning and disinfection

To clean and disinfect the attachments either of the two following methods may be used:

Manual Method

- Place the attachments in a perforated tray or basket.
- Prepare an ultrasonic bath and a detergent and disinfectant solution in accordance with the manufacturer’s specifications of the product used. Nueva Galimplant S.L.U. has validated this process using the product neodisher® septo PreClean diluted at 1% to prepare the ultrasonic bath.
- Immerse the basket with the attachments in the ultrasonic bath and keep them at a temperature of 20° C.
- Remove the attachments from the ultrasonic bath and rinse with plenty of purified water.
- If necessary, repeat the cleaning and disinfection step.

Automatic Method

- Place the attachments in the thermal washer-disinfector’s tray or basket. Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.
- Run the cleaning, disinfecting and drying cycle. Nueva Galimplant S.L.U. has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:
 - Cold water soak for 4 minutes.
 - 55° C clean for 10 minutes with detergent.
 - Cold water neutralisation for 6 minutes.
 - Cold water mid-cycle rinse for 3 minutes.
 - 93° C thermal disinfection for 5 minutes.
 - 110° C dry for 25 minutes.

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To guarantee optimum protection of the attachments we recommend that purified water is used throughout the whole process.

Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfector manufacturers.

Packaging

Once they are clean, disinfected and completely dry, the attachments must be packaged in a sterilisation bag.

The packaging of thermally sterilised medical devices must comply with the following requisites:

- Be appropriate for steam sterilisation, that is, have a minimum thermic resistance of 137° C and be sufficiently permeable for steam.
- Protect the devices sufficiently from mechanical damage.
- Comply with the EN ISO 11607 normative.

Sterilisation

- Load the autoclave with the conditioned attachments in the sterilisation bags. Do not overload the autoclave and always follow the equipment manufacturer’s instructions.
- Just like in the other steps of these instructions, we recommend the use of purified water to guarantee optimum quality.
- Nueva Galimplant S.L.U. recommends using an autoclave that complies with the requisites of the EN ISO 17665 normative, in addition to being validated and having adequate maintenance. The validation of the sterilisation was carried out in a steam autoclave using a 134° C pre-vacuum cycle for 6 minutes.
- To control the sterilisation process, the use of chemical indicators is recommended, such as, for example, the verification strips which the sterilisation bags come equipped with and which change colour when the sterilisation process has finished.

Note: it is important to always follow the autoclave manufacturer’s instructions and keep maintenance up to date.

Caution! The re-use and/or re-sterilisation of prosthetic attachments can cause risk of infection and/or cross contamination.

WARNINGS

Only original impression copings should be placed in the galimplant® implants.

The prosthetic and clinical practitioner must verify that the materials and techniques to be used during prosthetic and clinical procedures are appropriate and compatible with the galimplant® impression copings and must follow the manufacturer’s instructions for use of these materials.

Before the prosthetic part, the physician should carry out an exhaustive patient study, in order to determine a correct planning of the case.

For the placement of the galimplant® impressions coping, use only original wrenches of the galimplant® brand.

There is a risk of aspiration or ingestion of the products when they are used, for this reason, professionals should only use them.

The prosthetic abutments must solely and exclusively be placed during the implants osteointegration or if the implants have a great primary stability. Failure to observe the procedures set out in these instructions could harm the patient and give rise to one or all of the following complications:

- Aspiration or swallowing of the components.
- Infection.
- Damage to the implant, abutment, components or instruments.
- Implant failure.

Inform Nueva Galimplant of any complication that arises in relation to its products.

The practitioner must provided the patient with the information about side effects, interactions and complications in relation to galimplant® impressions coping.

Procedures for the disposal of the product

Once they have been extracted from the oral cavity, the prosthetic components must be disposed of as biological waste, in accordance with current local laws.

Additional information

For more information concerning the use of galimplant® products, contact your local representative or visit our website www.galimplant.com.

Validity

These instructions for use will be valid as May 2022.

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