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Marcado de Conformidad CE (ON 0051).

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Modelo

**IMPORTANT! Read the information contained in these operating instructions carefully before using the product.**

**NUEVA GALIMPLANT** complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical products. The galimplant® system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the galimplant® system is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

**Safety, responsibility and guarantee**

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U. The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

**Assistance and training**

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations. Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their obligation to personally make sure that the product is appropriate for the purposes and procedures planned. Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about its products through the platform: galimplant.training.

**Precautions and warnings**

Galimplant® drill system have not been sterilised before delivery, so they must be cleaned, disinfected and sterilised before use and it is essential that they are handled in an aseptic manner. The following descriptions are insufficient for the immediate application of any galimplant® surgical drills. This system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of the product in order to use them safely and correctly following the indications of these instructions for use.

**Patient Groups and Clinical Benefits**

The galimplant® surgical drills for guided surgery are intended for use in totally or partially edentulous people or with poor prognosis or non-restorable teeth that request treatment and that meet the criteria for psychological, medical, anatomical, surgical, prosthodontic and maintenance selection specified in these instructions for use. The clinical benefit of treatment is the replacement of the patient's missing teeth.

**Intended use**

The galimplant® surgical drills for guided surgery are designed for preparing the implant bed for the sole and exclusive fitting of a galimplant® implant using a surgical guide designed by galimplant3D.

**Indications**

The galimplant® surgical drills for guided surgery are indicated for the insertion of galimplant® dental implants through guided surgery. The use of components manufactured by third parties may cause mechanical and/or instrumental faults, tissue damage or unsatisfactory results.

**Description**

The galimplant® surgical drills system for guided surgery is manufactured using surgical steel in accordance with ISO 7153-1 and it is designed to be used together with rotational instruments with Type 1 connection according to ISO 1797. The galimplant® surgical drills system for guided surgery incorporate a top to ensure that the drilling depth does not exceed the length of the galimplant® implant.

In order to achieve the optimum primary stability of the implant, the following drilling sequence table must be adhered to. Initially it is used the pin drill to placing anchor pins of the surgical guide, then using the Galimplant® surgical drills range for guided surgery with the same length as the implant that will be insert, continue with a drilling sequence of increasing diameter as shown below:

Pin drill	Ø Drill (mm)			
	2	2,6	3,2	3,6
3,2 mm	•	•		
3,5 mm	•	•	•	
4 mm	•	•	•	•

The doctor must verify that the drilling sequence is suitable according with patient bone, that could be modified.

**Contraindications**

The relevant contraindications are those relating to surgical implant procedures in patients who are medically unfit for an oral surgical procedure or who are allergic or hypersensitive to stainless steel.

**Side Effects and complications**

The applicable side effects are those relating to surgical procedures to the placement of dental implants also those of the product that are described below:  
 Transient effects: localized infection, heat generation, overheating, burns, thermal damage.  
 Long-lasting effects: osteointegration problems, bone necrosis.  
 Complications: failure of the drill, lack of traceability, deformation, loss of surface, wear and loss cutting capacity.

**Storage and transport**

The galimplant® surgical drills are delivered packed in peel-pack bags inside a closed cardboard blister on which the label is placed and includes all product information, including the date of manufacture, which can be deduced from the batch number, the first four digits corresponding to the year and the next two to the month of manufacture (YYMMXX). This packaging contains and protects the product during storage and transport and does not have a specified shelf-life. Drills must be stored in a dry place in the original packaging at room temperature and not exposed to direct sunlight. When removing the product from its packaging one must verify that it is not damaged and that the product fits the description on the label.

**Surgical procedure**

**Planning**

Before surgery, an exhaustive patient evaluation, diagnostic and treatment planning process must be carried out, being the exclusive responsibility of the product user. Errors in treatment planning can cause the loss of implants. Prior to surgery, make sure the correct condition of the drills to use (cutting capacity and cleaning conditions) and that the hand-pieces, turbines and counter angles, where the drill will be connected are in perfect condition, both aseptic and functionally. Always check that the drill is correctly secured to the counter-angle. Prior to use on the patient, the correct turning of the drill should be verified. It is necessary to hold our products in order to help to avoid their aspiration during intraoral use. As well as following all of the obligatory precautions for all surgical procedures and guaranteeing the aseptic conditions, in this case, anatomical knowledge as well as assessment of the pre-operative radiographies should be used to ensure that no damage is caused to the nerves and veins damage during the drilling of the maxillary/mandibular bones. Take particular care of the patient's anatomical structures. The soft and hard tissues should be treated in a non-aggressive manner in order to create optimum conditions for the cicatrization of the implant. The implant bed should be prepared with extreme care. A heat trauma will hinder or prevent the healing of the dental implant. For this reason the following measures should be applied in order to reduce as far as possible the prevalence of excessive temperatures:

- Use the appropriate drilling speed for each drill diameter. A speed of between 300 and 800 rpm is recommended.

- Use sharpened instruments, each piece should be inspected following each use and cleaning, and the cutting edges should be checked, the drill should be disposed of if there is any evidence of wearing. The drills should be replaced after each 20 usages or when the cutting capacity has reduced.
  - An intermittent boring and drilling technique should be applied.
  - The drills should be highly refrigerated using pre-cooled (5°) sterile physiological saline solution (NaCl).
  - Use drills with increasing diameter (Ø 2 / 2.6 / 3.2 / 3.6 / 4.2 / 4.7 mm). We recommend that the diameter of the final drill should be 0.4mm less than the diameter of the implant, this value may vary depending on the type of bone.
- NOTE:** it is recommended to keep a spare set of drills in the event replacement is needed during a surgery.

**Cleaning, disinfection and sterilisation**

The galimplant® surgical drills are reusable and are sold in a NON-STERILE state. In accordance with the EN ISO 17664 Normative, Nueva Galimplant will provide all of the information necessary to process their health products labelled as "Non sterile" and which need to be cleaned, disinfected and sterilized prior to use, both for the first use and for all subsequent usages. All models of the galimplant® surgical drill are included in said health products.

**Warnings:**

- Do not reuse the drills without prior cleaning, disinfection and sterilisation. The new "non-sterile" products should be processed in the same way as the used products.
- For their own safety, staff should use adequate protective clothing and equipment at all times: gloves, mask, protective glasses, etc;
- Nueva Galimplant recommends the use of aldehyde free cleaner and disinfectants with proven effectiveness, in order to avoid the fixation of blood proteins (included on the VAH/DGHH list or which have the CE Marking).
- The use of neutral or slightly alkaline cleaning agents is recommended by Nueva Galimplant.
- It is essential to verify that the cleaning product used is compatible with the medical devices are going to processed and respect the specifications of use (temperature, concentration and time) indicated by the manufacturer.
- Nueva Galimplant cannot control the clinical manipulation procedures, the cleaning methods or any other user level conditions, therefore they will not be held liable for the sterilization of the products even when the recommended guidelines contained in this instruction manual have been followed.

**Limitations and restrictions:**

The repetition of the indicated processes has a minimal effect on the galimplant® surgical drills. The useful life of these is normally determined by the wear and tear and damage caused due to repeated use. The best way to assess the state of deterioration is by carrying out a thorough and functional check prior to use. Worn down cutting edges or signs of corrosion are indications that said drill is not suitable for use and it should be disposed of. The surgical drills should not be resharpened.

**Preparation in the place of use**

As soon as possible following their use, the galimplant® surgical drills should be transferred to the place where they are to be cleaned, in order avoid the drying-on of any blood residues, tissues, secretions or bone fragments. If it is likely that the cleaning process will be delayed, these should be covered with a moist cloth or stored in an airtight container filled with water in order to delay the drying-on of the organic material and to avoid the risk of contamination. To avoid contamination, the contaminated drills should not be put back into the surgical box. The used devices should be placed in a container that has been set up for the pre-cleaning which shall be carried out as per the following indications.

**Preparation before cleaning**

This stage must be carried out regardless as to whether the cleaning and disinfection processes are carried out manually or automatically.

- As soon as possible following their use, immerse the galimplant® surgical drills in a container that contains a cleansing solution with proteolytic enzymes, prepared in accordance with the manufacturer's specifications, respecting the temperature, concentration and immersion time. These enzyme solutions separate proteins and ensure that blood and other secretions do not dry on the devices.
- Make sure that all the drill's surfaces are covered by the solution and that the instruments do not touch each other.
- Manually brush with a nylon bristle brush (never metal, wool or steel brushes) until all of the visible dirt and residue has been eliminated, paying special attention to the zones which may not be reached by the brushing action, such as grooves or internal parts.
- Remove the drills from the container containing solution and rinse them thoroughly with plenty of water, preferably purified water, to ensure the reduction of organic material from mechanical drag.

**Cleaning and disinfection**

To clean and disinfect the drills either of the two following methods may be used:

**Manual Method**

- Place the galimplant® surgical drills in a perforated tray or basket.
- Prepare an ultrasonic bath and a detergent and disinfectant solution in accordance with the manufacturer's specifications of the product used.
- Nueva Galimplant has validated this process using the product neodisher® septo PreClean diluted at 1% to prepare the ultrasonic bath.
- Immerse the basket with the galimplant® surgical drills in the ultrasonic bath and keep them at a temperature of 20° C for 15 minutes.
- Remove the drills from the ultrasonic bath and rinse with plenty of purified water.
- If necessary, repeat the cleaning and disinfection step.
- Dry the devices using filtered compressed air and single use lint free cloth.

**Automatic Method**

Place the galimplant® surgical drills in the thermal washer-disinfector's tray or basket. Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.

- Run the cleaning, disinfecting and drying cycle. Nueva Galimplant has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:
  - Cold water soak for 4 minutes.
  - 55° C clean for 10 minutes with detergent.
  - Cold water neutralisation for 6 minutes.
  - Cold water mid-cycle rinse for 3 minutes.
  - 93° C thermal disinfection for 5 minutes.
  - 110° C dry for 25 minutes.

To guarantee optimum protection of the attachments we recommend that purified water is used throughout the whole process.

*Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfector manufacturers.*

**Inspection, maintenance and testing**

- Carefully inspect all of the galimplant® surgical drills in order to check if any damage has been caused to them (corrosion, loss of shape, breaking) or contamination. We recommend that the inspection is carried out in a well-lit environment or using a magnifying lens in order to improve visibility.
- If there are any traces of contamination on any of the devices the cleaning and disinfection process should be repeated.
- If there are any damage (deformation, failure of the drill, loss cutting capacity...), the drill should be discarded.

**Packaging**

Once they are clean, disinfected and completely dry, the drills must be packaged in a surgical box and then the box should be packaged in a sterilisation bag. The packaging of thermally sterilised medical devices must comply with the following requisites: - Be appropriate for steam sterilisation, that is, have a minimum thermic resistance of 137° C and be sufficiently permeable for steam. - Protect the individual devices or surgical boxes sufficiently from mechanical damage. - Comply with the EN ISO 11607 normative.

Once the devices have been packaged a label should be added to each sterilisation bag which contains data such as the sterilisation date and it should include the reference, batch and the number of devices if the bag does not have a transparent window which allows the user to see what is inside.

**Sterilisation**

- Load the autoclave with the conditioned galimplant® surgical drills in the sterilisation bags. Do not overload the autoclave and always follow the equipment manufacturer's instructions.
- Just like in the other steps of these instructions, we recommend the use of purified water to guarantee optimum quality.
- Nueva Galimplant recommends using an autoclave which complies with the requisites of the EN ISO 17665 normative, in addition to being validated and having adequate maintenance. The validation of the sterilisation was carried out in a steam autoclave using a 134° C pre-vacuum cycle for 6 minutes.

- To control the sterilisation process, the use of chemical indicators is recommended, such as, for example, the verification strips which the sterilisation bags come equipped with and which change colour when the sterilisation process has finished.

*Note: it is important to always follow the autoclave manufacturer's instructions and keep maintenance up to date.*

**WARNINGS**

Use the surgical drills only for the insertion of original galimplant® dental implants. The practitioner must provided the patient with the information about side effects, interactions and complications in relation to galimplant® surgical drills. The clinical practitioner must verify that the materials and techniques to be used during surgical procedure are appropriate and compatible with the galimplant® surgical drills and must follow the manufacturer's instructions for use of these materials. Failure to observe the procedures set out in these instructions could harm the patient and give

- rise to one or all of the following complications:
- Aspiration or swallowing of the components.
  - Infection.
  - Damage to the implant or instruments.
  - Implant failure.

In the event of a serious incident related to the use or result of use of the product, inform Nueva Galimplant (www.galimplant.com) and your National Competent Authority.

**Procedures for the disposal of the product**

If used, dispose the surgical drills as biological waste, according to the local regulations.


**Validity**

These instructions for use will be valid as of June 2024.

**BASIC UDI-DI**

843658232TF200903FIMPLSM

**SYMBOL**

-  Warning, read the instructions for use.
-  Non sterile product.
-  Do not use if packaging has been damaged.
-  Lot number.
-  Reference.
-  Consult the instructions for use.
-  Manufacturer.
-  CE (ON 0051) conformity marking.
-  Medical Device
-  Unique device identifies
-  Model

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