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**UDI-DI BÁSICO**

843658232TF200902PCICATH4

**IMPORTANT! Read the information contained in these operating instructions carefully before using the product.**

**NUEVA GALIMPLANT** complies with all the regulations established by the European laws and directives relating to the manufacture and distribution of medical products. The galimplant® implant system, attachments and instruments are certified and authorised for sale by a Notified Body, guaranteeing the perfect quality of all its products. The use and application of the galimplant® system is not under the manufacturer's control, and the user is responsible for any damages which may occur due to an incorrect manipulation or use of the product.

#### Safety, responsibility and guarantee

We guarantee the quality of our products within the framework of our sales and supply conditions. This product is part of a global concept and should only be used in conjunction with the associated original products according to the Nueva Galimplant instructions and recommendations. The non-recommended use of products manufactured by third parties together with galimplant® products will void any guarantee or other obligation, express or implicit, by Nueva Galimplant, S.L.U.

The user has sole responsibility for the procedure and use of the product. It is the practitioner's responsibility to use the product in accordance with the current instructions for use and to determine if it is appropriate for each patient's individual situation. In particular, this responsibility is expressly applicable in the case of its use in non-recommended procedures. If in doubt, the user should consult us.

#### Assistance and training

Technical assistance concerning the use of our products will be carried out verbally, by written means, by electronic means or by demonstrations. Said assistance will be given according to our knowledge of the current state of science and techniques when the product is on sale. This does not exempt the user from their obligation to personally make sure that the product is appropriate for the purposes and procedures planned.

Nueva Galimplant, with the aim of ensuring the correct use of its products, offers courses and practical sessions about its products through the platform: galimplant.training.

#### Precautions and warnings

In general terms, our products must be held to avoid their aspiration during intraoral use or insertion in the implant.

Galimplant® healing abutments have not been sterilised before delivery, so they must be cleaned, disinfected and sterilised before use and it is essential that they are handled in an aseptic manner.

Galimplant® healing abutments are of single use. The reusing of these types of products could lead to a deterioration of their characteristics which involves risk of tissue infection, surgical or prosthodontic failure and/or damage to the patient's health.

The following descriptions are insufficient for the immediate application of any abutment from the galimplant® system, this system should only be used by dentists, doctors and surgeons. The dentist/doctor/surgeon must have knowledge in oral implantology and experience in the manipulation of galimplant® products, in order to use them safely and correctly, following the indications of these instructions for use.

#### Patient Groups and Clinical Benefits

The galimplant® healing abutments with direct implant connection are intended for use in totally or partially edentulous people or with poor prognosis or non-restorable teeth that request treatment and that meet the criteria for psychological, medical, anatomical, surgical, prosthodontic and maintenance selection specified in these instructions for use. The clinical benefit of treatment is the replacement of the patient's missing teeth.

## Healing Abutments

#### Intended use

Our healing abutments with direct implant connection are intended for provisional use, solely and exclusively for the healing and configuration of soft tissues, to prevent the collapse of the tissue following surgical intervention.

#### Indications

The use of our healing abutments is recommended to prepare the emergence profiles of the mucous membranes before the prosthetic loading, thereby enabling the healing and configuration of the soft tissues around the implant's platform.

#### Description

Galimplant® healing abutments are made from grade 5 titanium and can be used in galimplant® implants with different connection types (internal and external), platforms, diameters and lengths.

See in the PRODUCTS section the implant models compatible with each of the configurations presented below.

Our galimplant® healing abutments are available for internal and external connections in two configurations:

- **Internal connection:**
  - Straight Healing Abutment, with heights of 1, 2, 4 and 6 mm.
  - Conical Healing Abutment, with heights of 4 and 6 mm.
- **External connection:**
  - Straight Healing Abutment, with heights of 1, 2, 4 and 6 mm.
  - Conical Healing Abutment, with heights of 4 and 6 mm.

The choice of the healing abutment height will depend on the height of the soft tissue.

Recommended torque 8-10 Ncm.

Do not re-use. Single use.

See the following table for the implant models compatible with each of the configurations.

| CONFIGURATIONS           |        |                     |                     |
|--------------------------|--------|---------------------|---------------------|
|                          | Height | External Connection | Internal Connection |
| Straight                 | 1mm    | PC 4010             | PCS 04010           |
|                          | 2mm    | PC 4020             | PCS 04020           |
|                          | 4mm    | PC 4040             | PCS 04040           |
|                          | 6mm    | PC 4060             | PCS 04060           |
| Conical                  | 4mm    | PCC 4040            | PCCS 04040          |
|                          | 6mm    | PCC 4060            | PCCS 04060          |
| Compatible Implant Model |        | ICE<br>IPXE         | ICI<br>IPX          |

## Healing Abutments

#### Intended use

The healing abutments for galimplant® multi-position abutments are intended for provisional use, solely and exclusively for the healing and configuration of soft tissues, to prevent the collapse of the tissue following surgical intervention.

#### Indications

The use of healing abutments for galimplant® multi-position abutments is recommended to prepare the emergence profiles of the mucous membranes before the prosthetic loading, thereby enabling the healing and configuration of the soft tissues around the implant's platform.

#### Description

The healing abutments for galimplant® multi-position abutments are made from grade 5 titanium and can be used in galimplant® multi-position abutments with different configuration types, connection (internal and external), heights and angulations.

Recommended torque 8-10 Ncm.

Do not re-use. Single use.

See the following table for the prosthetic abutment models compatible with each of the configurations.

| CONFIGURATIONS                             |                      |  |
|--|----------------------|--|
| COMPATIBLE ABUTMENT MODEL                  | External Connection  | Internal Connection                          |
| Multi-Position Aesthetic Straight Abutment | PCT 4030<br>PCT 4830 | PCT 4030<br>PCT 4830<br>PCT 5030<br>PCT 5830 |

|                                  |                      |                      |
|----------------------------------|----------------------|----------------------|
| Multi-Position Straight Abutment | PCM 4030<br>PCM 4830 | PCM 4030<br>PCM 4830 |
| Multi-Position Angled Abutment   | PCM 4030<br>PCM 4830 | PCM 4030<br>PCM 4830 |
| Compatible Implant Model         | ICE<br>IPXE          | ICI<br>IPX           |

#### Contraindications

The use of healing abutments is not recommended for patients:

- who do not meet the necessary medical requirements for dental implant surgery.
- who are allergic or hypersensitive to commercially pure titanium or Ti-6Al-4V titanium alloy (titanium, aluminium and vanadium).

#### Side effects, interactions and complications

The side effects, interactions and complications are those derived from implant placement procedures.

Possible complications following prosthetic restoration include:

- Complications: Loosening/fracture of the abutment and/or screw, mechanical problems associated with the prosthesis, loss of prosthesis retention, loss of the hole sealing of the screw-retained prosthesis, microgap in the implant-abutment connection and microleakage.

If one of these disorders/complications appears, the patient should consult their specialist doctor who will give them the corresponding instructions.

#### Principles for treatment planning

Before surgery, an exhaustive patient evaluation, diagnostic and treatment planning process must be carried out, being the exclusive responsibility of the product user. Errors in treatment planning can cause the loss of implants.

#### Storage and transport

The healing abutments are delivered packed in peel-pack bags inside a closed cardboard blister on which the label is placed and includes all product information, including the date of manufacture, which can be deduced from the batch number, the first four digits corresponding to the year and the next two to the month of manufacture (YYYYMMXX). This packaging contains and protects the product during storage and transport and does not have a specified shelf-life.

Healing abutments must be stored in a dry place in the original packaging at room temperature and not exposed to direct sunlight.

When removing the product from its packaging one must verify that it is not damaged and that the product fits the description on the label.

**Caution!** All our healing abutments are sold in a NON STERILE state.

Before use, they must be cleaned, disinfected and sterilised, observing the procedures validated by Nueva Galimplant SLU, included in the Cleaning/Disinfection/Sterilization section. Said procedures should be carried out before intraoral use. The repetition of the procedures described in this section do not alter the characteristics of these devices. The non-observance of these instructions may cause cross infections.

#### Cleaning/Sterilisation/Sterilization

The use of neutral or slightly alkaline cleaning agents is recommended. It is essential to verify that the cleaning product used is compatible with the medical devices are going to processed and respect the specifications of use (temperature, concentration and time) indicated by the manufacturer.

#### Cleaning and disinfection

To clean and disinfect the attachments either of the two following methods may be used:

##### Manual Method

- Place the attachments in a perforated tray or basket.
- Prepare an ultrasonic bath and a detergent and disinfectant solution in accordance with the manufacturer's specifications of the product used. *Nueva Galimplant S.L.U. has validated this process using the product neodisher® septo PreClean diluted at 1% to prepare the ultrasonic bath.*
- Immerse the basket with the attachments in the ultrasonic bath and keep them at a temperature of 20° C.
- Remove the attachments from the ultrasonic bath and rinse with plenty of purified water.
- If necessary, repeat the cleaning and disinfection step.

##### Automatic Method

- Place the attachments in the thermal washer-disinfector's tray or basket. Nueva Galimplant S.L.U. recommends thermal washer-disinfectors that comply with the EN ISO 15883 normative.
- Run the cleaning, disinfecting and drying cycle. *Nueva Galimplant S.L.U. has validated the Vario-TD program by the professional brand Miele, using the detergent neodisher® MediClean forte diluted at 0.5%. This program consists of:*
  - Cold water soak for 4 minutes.
  - 55° C clean for 10 minutes with detergent.
  - Cold water neutralisation for 6 minutes.
  - Cold water mid-cycle rinse for 3 minutes.
  - 93° C thermal disinfection for 5 minutes.
  - 110° C dry for 25 minutes.

To guarantee optimum protection of the attachments we recommend that purified water is used throughout the whole process.

Note: it is important to always follow the instructions provided by the detergent and thermal washer-disinfector manufacturers.

#### Packaging

Once they are clean, disinfected and completely dry, the attachments must be packaged in a sterilisation bag.

The packaging of thermally sterilised medical devices must comply with the following requisites:

- Be appropriate for steam sterilisation, that is, have a minimum thermic resistance of 137° C and be sufficiently permeable for steam.
- Protect the devices sufficiently from mechanical damage.
- Comply with the EN ISO 11607 normative.

#### Sterilisation

- Load the autoclave with the conditioned attachments in the sterilisation bags. Do not overload the autoclave and always follow the equipment manufacturer's instructions.
- Just like in the other steps of these instructions, we recommend the use of purified water to guarantee optimum quality.
- Nueva Galimplant S.L.U. recommends using an autoclave which complies with the requisites of the EN ISO 17665 normative, in addition to being validated and having adequate maintenance. The validation of the sterilisation was carried out in a steam autoclave using a 134° C pre-vacuum cycle for 6 minutes.
- To control the sterilisation process, the use of chemical indicators is recommended, such as, for example, the verification strips which the sterilisation bags come equipped with and which change colour when the sterilisation process has finished.

Note: it is important to always follow the autoclave manufacturer's instructions and keep maintenance up to date.

**Caution!** The re-use and/or re-sterilisation of healing abutments can cause risk of infection and/or cross contamination.

#### Documentation, traceability and implant card.

Nueva Galimplant S.L.U. recommends keeping an exhaustive clinical, radiological and photographic documentation. All products can be identified through their reference and batch code and / or UDI code. To maintain the mandatory traceability of these products, two adhesive labels are included, one for the medical record and the other for the patient, the latter must be attached to the implant card supplied with the products and which must be delivered to the patient fully completed once the surgery is over.

#### WARNINGS

Only original healing abutments should be placed in the galimplant® implants. The prosthetic and clinical practitioner must verify that the materials and techniques to be used during prosthetic and clinical procedures are appropriate and compatible with the galimplant® healing abutments and must follow the manufacturer's instructions for use of these materials.

Before the prosthetic part, the physician should carry out an exhaustive patient study, in order to determine a correct planning of the case.

For the placement of the healing abutments, use only original wrenches of the brand galimplant®. Consult the compatible product models in the catalogue available at [www.galimplant.com](http://www.galimplant.com).

There is a risk of aspiration or ingestion of the products when they are used, for this reason, they should only be used by professionals.

The prosthetic abutments must solely and exclusively be placed during the implants osteointegration or if the implants have a great primary stability.

Failure to observe the procedures set out in these instructions could harm the patient and give rise to one or all of the following complications:

- Aspiration or swallowing of the components.
- Infection.
- Damage to the implant, abutment, components or instruments.
- Implant failure.

To guarantee a good long term treatment result, the clinician must inform the patient about appropriate home oral hygiene and the importance of attending periodic clinical and radiographic consultations to perform a control and preventive diagnosis of any side effect and/or complication.

In the event of a serious incident related to the use or result of use of the product, inform Nueva Galimplant ([www.galimplant.com](http://www.galimplant.com)) and your National Competent Authority.

The practitioner must provide the patient with the information about side effects, interactions and complications in relation to galimplant® prosthetic abutments.

#### Procedures for the disposal of the product

Once they have been extracted from the oral cavity, the prosthetic components must be disposed of as biological waste, in accordance with current local laws.


#### Summary of Safety and Clinical Performance

The Summary of Safety and Clinical Performance (SSCP) for the product is available in the European database EUDAMED. The document can be obtained from the website: <https://ec.europa.eu/tools/eudamed>.

#### Validity

These instructions for use will be valid as of February 2025.

#### SYMBOLS

|   |  |
|---|--|
|    | Warning, read the instructions for use   |
|    | Non sterile product                      |
|    | Do not use if packaging has been damaged |
|    | Single-use device                        |
|    | Lot number                               |
|    | Reference                                |
|    | Consult the instructions for use         |
|   | Manufacturer                             |
|  | CE (0051) Conformity Marking             |
|  | Medical Device                           |
|  | Unique Device Identification             |
|  | Model                                    |

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